



User manual

AGING
SKILLFULLY!



Contents



- 02** *About EM-BRACE*
- 03** *Instructions*
- 04** *Safety Information*
- 05** *About the Bracelet*
- 06** *Bracelet Placement*
- 07** *About the Control Station*
- 08** *Warranty and Recycling*




About EM-BRACE

EM-BRACE is a set of two devices designed to monitor the elderly: a bracelet (a) and a control station (b). This product is suitable for people living alone who still have a degree of mobility. It aims to measure the user vital signs as well as the room parameters. The data can be verified at any time by authorized people through the EM-BRACE mobile application.



a) Bracelet



b) Control Station

Instructions

01

Remove the two devices from the package.

02

Put on the bracelet and adjust the strap.

03

Place the control station on a flat surface.

04


Recharge the devices whenever needed.

Safety Information

Read all safety information before using the device to ensure proper use.

Follow the warnings below to prevent incidents, such as fires or explosions, injuries, or damage to the devices:

- Do not expose the device to physical impact or damage.
- Use only EM-BRACE-approved batteries, chargers, and cables specifically designed for your device.
- Prevent the multipurpose jack and battery terminals from coming into contact with foreign materials such as metal, liquids, or dust.
- Do not disassemble the product.
- If any part of the product is broken, stop using the device immediately. Use the device again only after it has been repaired at an EM-BRACE Service Centre.
- Do not allow children or animals to chew or suck the devices.
- The device can be used in locations with an ambient temperature of 0 °C to 35 °C. You can store the device at an ambient temperature of -20 °C to 50 °C. Using or storing the products outside of the recommended temperature ranges may damage the device or reduce the battery's lifespan.
- Comply with all safety warnings and regulations regarding electronic device usage.




About the bracelet

The bracelet contains two sensors for monitoring pulse and temperature. A gyroscope is included to detect falls. In case of emergency, press the red button. The red LED will light up if a recharge is needed. Battery autonomy is about 2 weeks.

To charge the bracelet:

1. Remove the bracelet from your wrist.
2. Plug the charging cable into the USB port on your computer or a USB wall charger.
3. Connect the bracelet to the charger.




Bracelet Placement

For all-day wear, the bracelet should usually rest a finger's width below your wrist bone and lay flat, as you would normally wear a watch.

For improved heart rate accuracy do not wear the bracelet too tight. A tight band restricts blood flow, potentially affecting the heart rate signal. It should also be slightly tighter (snug but not constricting) during exercise than during all-day wear.






About the Control Station

The control station contains a sensor that measures the temperature, pressure, and humidity of the room. All the parameters registered by the two devices will be displayed on the screen. The data is stored, encrypted, and can be verified on the EM-BRACE mobile application.

To charge the control station:

1. Plug the charger into the outlet.
2. Connect the control station to the charger.




Warranty and Recycling

- This product has a 5-year warranty. Contact an authorized service center for maintenance.
- Please contact our customer service to recycle the product at the end of its use.
- Do not dispose of the product with other household waste.

Warranty information and the Return Policy can be found online at
www.em-brace-isep.com





Declaration of Conformity

EM-BRACE declares that this product complies with the essential requirements and relevant European Directives.

Copyright©2020 EM-BRACE